

Integrating Signaling Kinetics into GPCR Compound Profiling

Leigh A Stoddart¹, Julija Sirina¹, Nicola C Dijon², Nicholas D Holliday¹

¹ Excellerate Bioscience, Biocity, Nottingham UK

² Cell Signalling Research Group, School of Life Sciences, University of Nottingham, Nottingham, UK

- ❖ Optimisation of ligand binding and signalling kinetics is useful to consider in early stage compound profiling
- ❖ For G protein-coupled receptors (GPCRs) there is a relationship between the kinetics of signalling pathway activation and desired functional phenotype
- ❖ Common GPCR assay are endpoint based with responses measured at different timepoint
- ❖ Can new technologies be applied to develop robust kinetic signalling assay suitable for compound profiling?


❖ Excellerate Bioscience is a contract research organisation (CRO) specialising in molecular & cellular pharmacology and immunology.

❖ Based in BioCity Nottingham, our team of drug discovery scientists work on projects with leading pharmaceutical biotechnology & academic institutions around the globe.


❖ We combines the best of academic and industrial science to develop innovative approaches that will improve the efficiency and translatability of *in vitro* pharmacological profiling.

❖ Investing in capabilities, R&D methodology and assay platform development


NanoBiT® technology



cADDIS biosensor technology



- ❖ Circularly permuted GFP positioned between the catalytic and enzymatic domains of EPAC



- ❖ Luciferase complementation assay
- ❖ Luminescence measured in real time, in live cells at 37°C
- ❖ Cells stably expressing β_2 AR-LgBiT and β arrestin2-SmBiT





Figure 1 Potency and efficacy of reference β AR agonists in β arrestin2 recruitment assay

(a) Kinetic profile of agonist mediated β arrestin2-SmBiT recruitment to β_2 AR-LgBiT resulting in an increase in luminescence. Concentration response curves for (b) isoprenaline, (c) salbutamol or (d) BRL 37344 at 5 min (black circles) or 30 min (red squares). Concentration response curves can be plotted from any time point in the kinetic read. (e) Average potency values from 4 independent β arrestin2 recruitment assays.

	5 min		30 min	
	pEC ₅₀	Rmax (% 10 μ M iso)	pEC ₅₀	Rmax (% 10 μ M iso)
Isoprenaline	7.29 ± 0.25	100	7.77 ± 0.35	100
Salbutamol	6.50 ± 0.23	29.8 ± 10.4	6.89 ± 0.33	59.8 ± 4.9
BRL 37344	7.63 ± 0.17	4.53 ± 0.42	7.37 ± 0.20	14.2 ± 1.1



- ❖ Fluorescence measured in real time, in live cells at 37°C
- ❖ cADDIS sensor contained in BacMam vector and transduced into cells stably over-expressing β_2 AR




Figure 3 Potency and efficacy of reference β AR agonists in β arrestin2 recruitment assay

(a) Kinetic profile of cAMP levels in cells expressing SNAP- β_2 and cADDIS downwards sensor upon β_2 stimulation with agonist. Concentration response curves for (b) isoprenaline, (c) salbutamol or (d) BRL 37344 at 3 min (black circles) or 30 min (red squares). Concentration response curves can be plotted from any time point in the kinetic read. (e) Average potency values from 2 independent experiments.

	5 min	30 min
	pEC ₅₀	pEC ₅₀
Isoprenaline	8.03 ± 0.20	9.18 ± 0.44
Salbutamol	7.52 ± 0.11	8.49 ± 0.02
BRL 37344	6.61 ± 0.59	7.63 ± 0.04




Figure 2 Miniaturization of β arrestin recruitment assay

Cells expressing β_2 AR-LgBiT and β arrestin2-SmBiT were (a) grown in an adherent monolayer in a 96-well plate or (b) grown in a T75 culture flask, removed from the flask non-enzymatically and added in suspension to a 384-well plate. Cells were incubated with furimazine prior to the addition of increasing concentrations of 4 reference agonists. Data shown is 5 min after the addition of agonist and is representative of one experiment performed in duplicate.




Figure 4 Long term cAMP kinetics

(a) HEK293 cells were transduced with cADDIS downwards sensor and 24h post-transduction stimulated with (red squares) and without (black circles) 10 μ M forskolin (FSK). Fluorescence was measured every 60 min for 8 h. After 8h, in the presence of FSK, the fluorescence increases indicating a desensitization of FSK-stimulated cAMP production. (b) To confirm responsiveness of the cells to FSK, cells were treated with 10 μ M FSK after 8h or 24h in assay buffer.